

## **Appendix 2: AutoGen.Net Instructions** - Firmware Download Procedure

Patent # 5,841,675


New firmware releases for the OES LPA can be downloaded into the unit using the following procedure. The current firmware version can be determined when viewing WIRESKAN.NET Software.

To download firmware to the monitor, these materials are needed:

1. A laptop, PDA or PC with the OES "AutoGen.NET" program loaded.
2. The updated firmware file (filename has the form: \*.S)
3. A serial communication cable (OES part #: AU008A)

### **Procedure for loading the software into the Laptop or PC:**


Power on your PC, once the PC is powered on insert CD provided with the AutoGen.NET version install on it. Run windows explorer program and select the CD drive, then select "Autogen.NET Installer".



Double click on this icon and the Autogen software will be installed onto your PC. An "OES, Inc" directory is automatically created. Select Autogen.Net from the directory and the following window will appear.

Open the file to be downloaded. (The file will be an "S record file").

Confirm the serial connection and com port between your PC and the LPA.



Click on "START" and the file will be downloaded. A successful download is confirmed. If the download is unsuccessful, re-check settings and download again.

